



Quick Start Guide

Mirus-Connect: Tubing Set for Mirus Evo™ Nanopump

Mirus-Connect: Tubing Set for Mirus Evo Nanopump

Mirus connect tubing set for the Mirus Evo nanopump contains:

- 1 pump inlet cable
- 1 pump outlet cable
- 1 pump-to-MultiFlow8 cable

To connect the tubing to the Mirus Evo nanopump, follow the steps illustrated below.



Connect inlet cable from flask to pump as shown



Connect outlet cable as shown



Connect cable to MultiFlow8 as shown

